



The closest you will come to "Plaster" in a Synthetic Cast Tape

- Conforms nicely to complex shapes
- Allows for the practitioner to produce a more conforming cast
- No plaster mess
- Ample time to work with before it hardens
- No need for special glove, a regular casting glove is fine to use
- Resin does not stick to the glove like typical fiber glass casting tape
- Very rigid cast once finished
- Cost effective (based on comparable conforming casting tapes)



Product Description	Product #
C-Form Casting Tape - White 2" x 4yd (box=10 rolls)	4013-2, 10
C-Form Casting Tape - White 3" x 4yd (box=10 rolls)	4013-3, 10
C-Form Casting Tape - White 4" x 4yd (box=10 rolls)	4013-4, 10
C-Form Casting Tape - White 5" x 4yd (box=10 rolls)	4013-5, 10

C-Form Usage Instructions:

Needed items:

1. Container of water.
2. Cutting strip or cutting channel for orthotic casting, or unusual shapes
3. EZ Cutter, or scissors
4. 2 nylon hose, or 2 stockinette
5. Gloves
6. Towel

Tip #1: The first layer of C-FORM Tape should be rolled gently. Then gently stretch second layer to release resin from casting tape.

Tip #2: When you feel the heat from the resin setting off, this is when you want to start to shape the cast to fit the patient.

Tip #3: Use HOT water for faster activation.

Tip #4: Use THE EZ CUTTER or a sharp pair of scissors when removing the cast.