

# 1271S / 1271T

## Double Pyramid Adapter, Ti

The 1271S or 1271T offers a double pyramid on both ends, available in titanium or stainless steel, and provides excellent adjustability options especially with going from receiver to receiver with limited space.

- Weight: 139g (SS) / 77g (Ti)
- Build Height: 10mm / 0.4in
- Weight Limit: 136kg / 300lbs (SS)  
170kg / 375lbs (Ti)
- Ideal for low clearance situations requiring a pyramid receiver



### 1. Application

The 1271, Double Pyramid Adapter, is to be exclusively used for the prosthetic fitting of lower extremity amputations.

### 2. Description

The 1271 is used to connect any pyramid receiver to another pyramid receiver. The 1271 has an overall 10mm / 0.4in build height. position.

Apply Loctite to all screws once the alignment and connection is finalized.

### 3. Warranty

Warranty applies only when the product is used as intended without modification, following all manufacturer's recommendations. ST&G offers 2 years warranty.

### 4. Environment:

Store at room temperature, and avoid extremely hot, cold, or moist storage environments Storage Temperature Range: between temperatures of -20°C to 80°C [-4°F and 176°F]

Do not store in wet or very hot environment.



### 5. CE Conformity

This product meets the requirements Council EU 2017/1745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future use.



### \*Notice\*

Cross tighten the set screws, which are side by side, using a torque wrench set to 15 Nm.