



Pediatric SACH Foot Adapter P231A

Instructions For Use



P231A

1. Application & Description

- This adapter is exclusively for lower extremity prosthetic fittings.
- It is a pediatric 22mm pyramid receiver system to a pediatric SACH foot utilizing a Foot Bolt.
- The four set screws anchor the mating system and allow for angular alignment.

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
P231A	Aluminum	7mm	28g	55kg/122lbs

- **System Size: 22mm**
- **Recommended Torque:** Always cross tighten Alignment Set Screws - loctite all fasteners
- **M6 Alignment Set Screw Torque:**
 - (Aluminum): **6Nm**
- **M8 Cap Screw Foot Bolt Torque:**
 - (Aluminum): **20Nm**

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F)
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet; very hot or extremely cold environments.
- **Warranty:** 1 years, provided it is used as intended without modification and only the original user.
- Do not reuse on another prosthesis
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



P231A IFU REV B (9-30-2026)

Specifications subject to change.