



# Pediatric Tube Clamp P212A

## Instructions For Use



P212A

### 1. Application & Description

- This adapter is exclusively for lower extremity prosthetic fittings.
- It connects any pediatric 22mm pyramid system to a 22mm pylon.
- The four set screws anchor the mating system and allow for angular alignment.

### 2. Technical Specifications

| P/N   | Material | Build Height | Weight | Weight Limit |
|-------|----------|--------------|--------|--------------|
| P212A | Aluminum | 13.5mm       | 31g    | 55kg/122lbs  |

- **System Size: 22mm**
- **Recommended Torque:** Always cross tighten Alignment Set Screws
- **M6 Alignment Set Screw Torque:**
  - (Aluminum): **6Nm**
- **M5 Cap Screw Torque:**
  - (Aluminum): **5Nm**
- **Thread Locker:** Apply Loctite once alignment and connection are finalized

### 3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F) 
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet; very hot or extremely cold environments.
- **Warranty:** 1 years, provided it is used as intended without modification and only the original user.
- Do not reuse on another prosthesis 
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.

P212A IFU REV B (9-30-2026)

Specifications subject to change.