



Pediatric Tube Clamp Assembly 22mm, Al

Instructions For Use

P211A

P211A



1. Application & Description

- This 22mm adapter is exclusively for Pediatric lower extremity prosthetic fittings.
- It connects any pediatric pyramid system to a tube clamp, with four set screws that allow for angular alignment.

2. Technical Specifications

P/N	Material	Pylon Length	Weight	Weight Limit
P211A	Aluminum	270mm	211g	55kg/122lbs

- **Compatibility:**
 - The 270mm Pediatric Tube Clamp Assembly, Al is designed for use only ST&G 22mm tube clamps.
- **Installation Requirements:**
- **Recommended Torque:** (Always cross tighten Alignment Set Screws)
- **Alignment Set Screw: 6Nm**
- **Thread Locker:** Apply Loctite once alignment and connection are finalized.

CAUTION: Do not clamp the tube in a vise or similar device, as this may damage the tube

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F)
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 1 years, provided it is used as intended without modification, and only for the original intended user
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



P211A IFU REV B (9-30-2026)

Specifications subject to change.