



Foam Block, 4-Hole Adapter

1209W

Instructions For Use



1209W

1. Application & Description

- The 1209W, Foam Block 4-Hole Adapter is to be exclusively used for the prosthetic fitting of lower extremity amputations.
- It is suitable for both Transtibial and Transfemoral sockets connecting a 4-Hole adapter to the socket.

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
1209W	Foam	69mm	290g	136kg/300lbs

- **M6 Flat/ Button Head Screw Torque: 10Nm**

3. Installation Overview

- **Attachment:** The block is designed to be laminated into a socket, and typically in a two-step lamination. Socket is bench aligned and bonded onto the block, (the block surface must be sealed with resin first) with a very strong bonding agent. Once alignment is set, trim and blend the block for second lamination.
- **Compatibility:** Connects any modular component featuring a 4-hole system to a socket
- **Fastening:** Tighten all screws that are side by side using a torque wrench

4. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F).
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



1209W IFU REV. E (3-26-2026)
Specifications subject to change.