



EZ Adapt Adjustable Height Tube Clamp Pyramid Receiver, Heavy Duty, Ti

1753T

Instructions For Use



1753T

1. Application & Description

- The Variable Height Double Ended Adapter is to be exclusively used for the prosthetic fitting of lower extremity amputations to connect any pyramid system to a pylon.
- The 1713T provides:
 - **12.7mm / 0.5in of height adjustment**
- **Compatible with 30mm pylons**
- **Description:** This titanium adapter features a female pyramid receiver with a 12.7mm offset threaded adapter designed for high strength and versatility when space is limited.

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
1753T	Titanium	45mm to 57mm	203gm	125kg/275lbs

- **Recommended Torque:** Always cross tighten Alignment Set Screws
- **M8 Alignment Set Screw Torque:**
 - (Titanium): **15Nm**
- **M6 Cap Screw Torque:**
 - (Titanium): **7Nm**

3. Installation & Maintenance

- **Alignment:** After alignment is complete, set screws must be torqued to **15Nm**.
- **Thread Locker:** Apply Loctite or a similar thread locker to all screws during final assembly
- **Cross Torque:** Use a torque wrench to cross-torque all screws to ensure even distribution

4. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F) 
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user
- Do not reuse on another prosthesis. 
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.

1753T IFU REV. G (3-20-2026)

Specifications subject to change.