



Double Ended Adapter 1201A, 1202A, 1203A, & 1204A

Instructions For Use



1201A



1202A



1203A



1204A

1. Application & Description

- This adapter is for use in extremity prosthetic fittings.
- **Strictly for K2 only.**
- It connects any pyramid system on both ends of the receiver.
- The four set screws anchor the mating system and allow for angular alignment.

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
1201A	Aluminum	32mm	100g	158kg/350lbs
1202A	Aluminum	45mm	124g	158kg/350lbs
1203A	Aluminum	60mm	145g	158kg/350lbs
1204A	Aluminum	75mm	167g	158kg/350lbs

- **Recommended Torque:** Always cross tighten Alignment Set Screws
- **M8 Alignment Set Screw Torque:**
 - (Aluminum): **10Nm**
- **Thread Locker:** Apply Loctite once alignment and connection are finalized

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F).
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



1201A_1202A_1203A_1204A IFU REV. B (5-26-26)

Specifications subject to change.