



30mm Carbon Pylon

Instructions For Use

1209C



1209C

1. Application & Description

- This **30mm** pylon is exclusively for lower extremity prosthetic fittings up to 375lb max weight limit for Trans-Tibial fittings. **Only compatible with ST&G 3212T Tube Clamp.**
- The included tubular inserts are to be used for clamping within the 3212T Tube Clamp. Failure to include this can lead to potential failure to sustain clamping forces.

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
1209C	Carbon Pylon	305mm	79g	170kg/375lbs

- **CAUTION:** Do not clamp the tube in a vise or similar device, as this may damage the tube

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F).
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



4. Reinforcement Recommendations

- The clinician is responsible for determining the adequate layup and reinforcement based on the specific patient's usage, activity level, and expected loading.

1209C IFU REV. F (3-26-2026)
Specifications subject to change.